


## Installation Instructions

### BED-BAR Model # 1200H

The BED-BAR can be attached to either the bed frame or the mattress frame. When attached to the bed frame the BED-BAR will remain stationary when the head of the bed is raised and lowered. If attached to the mattress frame the BED-BAR will move up and down when the head of the bed is raised and lowered.



1. Match the underside bracket with the corresponding bracket attached to the BED-BAR.
2. Rest the BED-BAR on the frame of the bed at desired attachment location.
3. Place bolts through the slots in the brackets (Note Figure 2) and tighten.